



ENV 62 | Polynuclear Aromatic Hydrocarbons (PAH) Analysis on BPX50

| | | | |
|-------------------------|------------------------------|-----------------------------|----------------------|
| Column Part No.: | 054751 | Flow: | On |
| Phase: | BPX50, 0.25 µm film | Average Linear Velocity: | 39 cm/sec at 65 °C |
| Column: | 30 m x 0.25 mm ID | Mode of Injection: | splitless |
| (PAH) standard: | 10 ng/ µL in dichloromethane | Purge on Time: | 0.5 min. |
| Initial Temp.: | 65 °C, 0.5 min | Purge on (split) Vent Flow: | 60 mL/min |
| Rate 1.: | 25 °C/min to 140 °C | Injection Volume: | 0.2 µL |
| Rate 2.: | 10 °C/min to 325 °C | Injection Temp.: | 250 °C |
| Final Temp.: | 325 °C, 15 min | Autosampler: | No |
| Detector Type: | MSD | Liner Type: | 4 mm ID Double Taper |
| Carrier Gas: | Helium, 9.7 psi | Liner Part Number: | 092018 |
| Carrier Gas Flow: | 1.1 mL/min constant | | |




ARO 08 | Analysis of Polynuclear Aromatic Hydrocarbons on HT8

| | | | |
|-------------------------|-------------------|--------------|---------------|
| Column Part No.: | 054462 | Rate: | 4 °C/min |
| Phase: | HT8, 0.25 µm film | Final Temp: | 380 °C, 5 min |
| Column: | 25 m x 0.22 mm ID | Carrier Gas: | He, 20 psi |
| Initial Temp: | 150 °C, 1 min | Detector: | FID |



GC Columns and Applications







GC Columns and Applications

AP 0040C | HT8: The Perfect PCB Column

Separation of CB31 & CB28



Chromatogram on the left clearly demonstrates the significant difference in selectivity of the HT8 column. By GC/MS, quantitation of CB28 using a standard 5% phenylpolysiloxane column is impossible as coelution with CB31 (with the same number of chlorines) occurs.



HT8 separates the two congeners by a full minute allowing quantitation to be performed with ease.

AROCLOR 1242

Column Part No.: 054676

Phase: HT8, 0.25 µm film

Column: 50 m x 0.22 mm ID

Initial Temp: 80 °C, 2 min

Rate 1: 30 °C/min

Temp 2: 170 °C


Rate 2: 3 °C/min

Final Temp: Split, 300 °C



Carrier Gas: He, 40 psi

Detector: ECD, 330 °C

| Congener # | Cl Position | Cl # | Identification by GC/MS |
|------------|-------------|------|-------------------------|
| 42 | 23-24 | 4 | ✓ |
| 96 | 236-26 | 5 | ✓ |
| 35 | 34-3 | 3 | ✓ |
| 64 | 235-4 | 4 | * |
| 72 | 25-35 | 4 | * |
| 103 | 246-25 | 5 | ✓ |
| 71 | 26-34 | 4 | ✓ |
| 41 | 234-2 | 4 | ✓ |
| 68 | 24-35 | 4 | ✓ |
| 37 | 34-4 | 3 | ✓ |
| 100 | 246-24 | 5 | ✓ |



| Congener # | Cl Position | Cl # | Identification by GC/MS |
|------------|-------------|------|-------------------------|
| 42 | 23-24 | 4 | ✓ |
| 96 | 236-26 | 5 | ✓ |
| 35 | 34-3 | 3 | ✓ |
| 64 | 235-4 | 4 | * |
| 72 | 25-35 | 4 | * |
| 103 | 246-25 | 5 | ✓ |
| 71 | 26-34 | 4 | ✓ |
| 41 | 234-2 | 4 | ✓ |
| 68 | 24-35 | 4 | ✓ |
| 37 | 34-4 | 3 | ✓ |
| 100 | 246-24 | 5 | ✓ |




TP-0187-C | Analysis of Polynuclear Aromatic Hydrocarbons on BPX50

| | | | |
|-------------------------|-----------------------------|----------------------|-------------------------------|
| Column Part No.: | 054701 | Temperature Profile: | Hold 70 °C for 1 min |
| Phase: | BPX50, 0.25 µm film | | 70 °C to 140 °C at 25 °C/min |
| Column: | 30 m x 25 µm ID | | 140 °C to 250 °C at 15 °C/min |
| Gas Flow: | 1.5 ml/min Helium | | 250 °C to 310 °C at 4 °C/min |
| Injection: | Split 1 µl (1 ng on column) | | Hold 310 °C for 8 min |
| Injection Temperature: | 250°C | | |




GC Columns and Applications



ENV 17 | Analysis of 16 Volatile Compounds in Drinking Water on BP624

| | | | |
|-------------------------|-------------------|-----------------|------------|
| Column Part No.: | 054826 | Final Temp.: | 170 °C |
| Phase: | BP624, 1.2 µm | Detector: | HP5870 MSD |
| Column: | 25 m x 0.22 mm ID | Injection Mode: | Splitless |
| Initial Temp.: | 50 °C, 2 min | Carrier Gas: | He, 15 psi |
| Rate: | 15 °C/min | | |






GC Columns and Applications

ENV 13 | Analysis of Volatiles from Drinking Water on BP624

| | | | |
|-------------------------|-------------------|-----------------|-------------------------|
| Column Part No.: | 054835 | | |
| Phase: | BP624, 3.0 µm | Rate 2: | 15 °C/min |
| Column: | 50 m x 0.53 mm ID | Final Temp.: | 210 °C, 1 min |
| Initial Temp.: | 35 °C, 2 min | Detector: | MSD, MJSC Jet Separator |
| Rate 1: | 8 °C/min | Injection Mode: | Purge & Trap |
| Temp 2: | 180 °C, 5 min | Carrier Gas: | He, 10 ml/min |

Note: Column which provides fast analysis of all EPA compounds. BP624 is also ideal for the analysis of many commonly used solvents.



Components

1. Carbon dioxide
2. Dichlorodifluoromethane
3. Chloromethane
4. Vinyl chloride
5. Acetaldehyde
6. Bromomethane
7. Chloroethane
8. Trichlorofluoromethane
9. Trichlorofluoroethane
10. Acrolein
11. Acetone
12. 1,1-Dichloroethene
13. Carbon disulfide
14. Methylene chloride
15. trans-1,2-Dichloroethene
16. Acrylonitrile

17. 1,1-Dichloroethane
18. Vinyl acetate
19. 2-Butanone (MEK)
20. cis-1,2-Dichloroethene
21. Bromochloromethane (Int. Std.)
22. 1,1,1-Trichloroethane
23. Carbon tetrachloride
24. 1,2-Dichloroethane-d4 (Surrogate)
25. 1,2-Dichloroethane
26. Trichloroethene
27. 1,2-Dichloroethene
28. Bromodichloromethane
29. 4-Methyl-2-pentanone
30. cis-1,3-Dichloropropene
31. Toluene-(d8) (Surrogate)
32. Toluene
33. trans-1,3-Dichloropropene
34. 2-Bromo-1-chloropropane (Int. Std.)

35. 1,1,2-Trichloroethane
36. 2-Hexanone
37. Tetrachloroethene
38. Dibromochloromethane
39. Chlorobenzene
40. Ethylbenzene
41. m,p-Xylene
42. o-Xylene
43. Styrene
44. Bromoform
45. 1,4-Dichlorobutane (Int. Std.)
46. Bromofluorobenzene
47. 1,1,2,2- Tetrachloroethene
48. 1,3-Dichlorobenzene
49. 1,4-Dichlorobenzene
50. 1,2-Dichlorobenzene




TP-0102-C | Analysis of Volatile Organic Pollutants on a Volatiles GC Column



GC Columns and Applications

| | | | |
|-------------------------|-------------------------------|--------------------------|----------------------|
| Column Part No.: | 054979 | Average Linear Velocity: | 35 cm/sec at 40 °C |
| Phase: | BPX-Volatiles 1µm film | Injection Mode: | Split |
| USEPA 502.2 mix: | 200 ppm in Methanol | Split Ratio: | 50:1 |
| Column: | 40m x 0.18mm ID | Injection Volume: | 1 µL |
| Initial Temp: | 40 °C, 0 min. | Injection Temperature: | 250 °C |
| Rate 1: | 6 °C to 210 °C | Autosampler: | No |
| Rate 2: | 15 °C to 250 °C | Liner Type: | 4 mm ID Single Taper |
| Final Temp: | 250 °C, 5 min | Liner Part Number: | 092017 |
| Detector Type: | Mass Spectrometer | Column Part Number: | 054979 |
| Carrier Gas: | He, 40.3 psi | ms-NoVent™ Part no.: | 113400 |
| Carrier Gas Flow: | 1.2 µL/min. | HP5973 restrictor: | 113409 |
| Constant Flow: | On | Full scan | 45-450 |



Notes. Chromatogram showing analysis of commonly screened volatile organic pollutants

Components

- | | | |
|-----------------------------|-------------------------------|---------------------------------|
| 1. Dichlorodifluoromethane | 20. Trichloroethene | 41. Bromobenzene |
| 2. Chloromethane | 21. 1,2-Dichloropropane | 42. 1,1,2,2-Tetrachloroethane |
| 3. Vinyl chloride | 22. Dibromomethane | 43. 1,2,3-Trichloropropane |
| 4. Bromomethane | 23. Bromodichloromethane | 44. n-Propyl benzene |
| 5. Chloroethane | 24. cis-1,3-Dichloropropene | 45. 2-Chlorotoluene |
| 6. Trichlorofluoromethane | 25. Toluene | 46. 1,3,5-Trimethylbenzene |
| 7. 1,1-Dichloroethene | 26. trans-1,3-Dichloropropene | 47. 4-Chlorotoluene |
| 8. Dichloromethane | 27. 1,1,2-Trichloroethane | 48. tert-Butylbenzene |
| 9. trans-1,2-Dichloroethene | 28. Tetrachloroethene | 49. 1,2,4-Trimethylbenzene |
| 10. 1,1-Dichloroethane | 29. 1,3-Dichloropropane | 50. sec-Butylbenzene |
| 11. 2,2-Dichloropropane | 30. Dibromochloromethane | 51. 1,3-Dichlorobenzene |
| 12. cis-1,2-Dichloroethene | 31. 1,2-Dibromoethane | 52. p-Isopropyltoluene |
| 13. Bromochloromethane | 32. Chlorobenzene | 53. 1,2-Dichlorobenzene |
| 14. Chloroform | 33. Ethylbenzene | 54. n-Butylbenzene |
| 15. 1,1,1-Trichloroethane | 34. 1,1,1,2-Tetrachloroethane | 55. 1,4-Dichlorobenzene |
| 16. 1,1-Dichloropropene | 35. p-Xylene | 56. 1,2-Dibromo-3-chloropropane |
| 17. Carbon tetrachloride | 36. m-Xylene | 57. 1,2,4-Trichlorobenzene |
| 18. Benzene | 37. o-Xylene | 58. Hexachlorobutadiene |
| 19. 1,2-Dichloroethane | 38. Styrene | 59. Naphthalene |
| | 39. Bromoform | 60. 1,2,3-Trichlorobenzene |
| | 40. Isopropylbenzene | |






GC Columns and Applications

ENV 20 | Analysis of Polychlorinated p-Dibenzodioxins on BPX5


| | |
|-------------------------|------------------------|
| Column Part No.: | 054114 |
| Phase: | BPX5, 0.25 μ m |
| Column: | 50 m x 0.22 mm ID |
| Initial Temp.: | 80 $^{\circ}$ C, 2 min |
| Rate 1: | 4 $^{\circ}$ C/min |
| Temp 2: | 220 $^{\circ}$ C |
| Rate 2: | 5 $^{\circ}$ C/min |

| | |
|-------------------|-----------------------------|
| Temp. 3: | 235 $^{\circ}$ C, 7 min |
| Rate 3: | 5 $^{\circ}$ C/min |
| Final Temp.: | 330 $^{\circ}$ C, 6 min |
| Detector: | High Resolution |
| Mass Spectrometer | He, 15 psi |
| Carrier Gas: | He, 300 psi |
| Injection Mode | Splitless, 270 $^{\circ}$ C |


Hexachlorodibenzodioxins




Tetrachlorodibenzodioxins




Heptachlorodibenzodioxins




Pentachlorodibenzodioxins



Octachlorodibenzodioxin



SGE wishes to acknowledge CARSO, 321 Avenue Jean Jaures, 69362 LYON CEDEX 7, FRANCE






GC Columns and Applications

ENV 57 | 8081 Organochlorine Pesticide Mix on BPX35

| | | | |
|-------------------------|------------------------------|-----------------------------|----------------------------|
| Column Part No.: | 054701 | | |
| Phase: | BPX35 0.25 µm film | Constant Flow: | On |
| Column: | 30 m x 0.25 mm ID | Average Linear Velocity: | 36 cm/sec at 40 °C |
| 8081 Standard: | 10 ng/ µL in dichloromethane | Injection Mode: | Splitless |
| Initial Temp.: | 40 °C, 1 min. | Purge on Time: | 1 min. |
| Rate 1: | 30 °C to 190 °C, 3 min | Purge on (Split) Vent Flow: | 60 mL/min. |
| Rate 2: | 10 °C to 300 °C | Injection Volume: | 1 µL |
| Final Temp.: | 300 °C, 5 min. | Injection Temp.: | 250 °C |
| Detector Type: | MSD | Autosampler: | No |
| Carrier Gas: | He, 10.0 psi | Liner Type: | 4 mm ID Double Taper Liner |
| Carrier Gas Flow: | 1.3 mL/min | Liner Part Number: | 092018 |



Components

1. 2,4,5,6-tetrachloro-meta-xylene
2. α-BHC
3. γ-BHC
4. β-BHC
5. Heptachlor
6. δ-BHC
7. Aldrin
8. Heptachlorepoxy
9. trans-Chlordane
10. cis-Chlordane
11. Endosulfan A
12. DDE
13. Dieldrin
14. Endrin
15. DDD
16. Endosulfan B
17. DDT
18. Endrin Aldehyde
19. Endosulfan Sulfate
20. Methoxychlor
21. Endrin Ketone
22. Decachlorobiphenyl

ENV 03 | Analysis of 18 Chlorinated Pesticides on BPX5


| | | | |
|------------------|-------------------|-----------------|---------------|
| Part No.: | 054125 | | |
| Phase: | BPX5, 0.5 µm film | Final Temp.: | 290 °C, 5 min |
| Column: | 25 m x 0.32 mm ID | Detector: | ECD at 310 °C |
| Initial Temp.: | 170 °C | Injection Mode: | Split |
| Rate: | 7 °C | Carrier Gas: | He, 7 psi |

Notes: Combined with the BPX608 column, BPX5 is the perfect confirmation column.

Components

20ng/ µL each component

1. α-BHC
2. γ-BHC
3. β-BHC
4. Heptachlor
5. δ-BHC
6. Aldrin
7. Heptachlorepoxy (isomer B)
8. Endosulfan I
9. 4,4'-DDE
10. Dieldrin
11. Endrin
12. 4,4'-DDD
13. Endosulfan II
14. 4,4'-DDT
15. Endrin aldehyde
16. Endosulfan sulfate
17. Methoxychlor
18. Endrin ketone







GC Columns and Applications

ENV 59 | 8141 Organophosphorous Pesticide Mix on BPX5

| | |
|-------------------------|------------------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| 8141 Standard: | 10 ng/ µL in dichloromethane |
| Initial Temp.: | 50 °C, 1 min |
| Rate 1: | 30 °C/min to 190 °C, 3 min |
| Rate 2: | 10 °C/min to 300 °C |
| Final Temp.: | 300 °C, 5 min. |
| Detector Type: | MSD |
| Carrier Gas: | He, 11.1 psi |
| Carrier Gas Flow: | 1.3 mL/min |

| | |
|-----------------------------|----------------------------|
| Constant Flow: | On |
| Average Linear Velocity: | 42 cm/sec at 50 °C |
| Injection Mode: | Splitless |
| Purge on Time: | 0.5 min |
| Purge on (Split) Vent Flow: | 60 mL/min |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Autosampler: | No |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |



Components

- | | | |
|-------------------------------------|----------------------|------------------------|
| 1. 4-Chloro-3-nitrobenzotrifluoride | 8. Naled | 17. Trichlorinate |
| 2. Dichlorvos | 9. Phorate | 18. Tetrachlorvinphos |
| 3. 1-Bromo-2-nitrobenzene | 10. Demeton | 19. Tokuthion |
| 4. α-Mevinphos | 11. Diazinon | 20. Impurity |
| 5. Tri-butylphosphate | 12. Disulfoton | 21. Fensulfothion |
| 6. Ethoprop | 13. Methyl parathion | 22. Impurity |
| 7. Sulfotepp | 14. Ronnel | 23. Triphenylphosphate |
| | 15. Chlorpyrifos | 24. Guthion |
| | 16. Fenthion | 25. Coumaphos |

ENV 45 | Organophosphorous Pesticides on BPX50

| | |
|-------------------------------|---------------------------------------|
| Column Part No.: | 054751 |
| Phase: | BPX50, 0.25 µm film |
| Mixture of: | 10 ng/ µL |
| Organophosphorous Pesticides: | 10 ng/ µL in |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp: | 50 °C , 1 min |
| Rate 1: | 30 °C/min to 200 °C, 3 min |
| Rate 2: | 10 °C/min to 310 °C |
| Final Temp: | 310 °C, 2 min |
| Detector Type: | FID, 320 °C |
| Carrier Gas: | He, 14.4 psi |
| Carrier Gas Flow: | 1.30 mL/min |
| Constant Flow: | On |
| Average Linear Velocity: | 30 cm/sec at 50 °C |
| Injection Mode: | Splitless |
| Purge On Time: | 0.5 min |
| Purge On (Split) Vent Flow: | 60 mL/min |
| Injection Volume: | 1.0 µL |
| Injection Temperature: | 240 °C |
| Autosampler: | Yes |
| Liner Type: | 4 mm ID FocusLiner™ with single taper |
| Liner Part Number: | 092003 |


| | |
|------------------------------|-----------------------|
| Column Part Number: | 054740 |
| Phase: | BPX50, 0.10 µm film |
| Mixture of 10 ng/ µL | 42 cm/sec at 50 °C |
| Organophosphorous Pesticides | Splitless |
| Column: | 10 m x 0.10 mm ID |
| Initial Temp.: | 70 °C , 1 min |
| Rate 1: | 25 °C/min to 320 °C |
| Rate 2: | N/A |
| Final Temp: | 320 °C, 0 min |
| Detector Type: | FID, 320 °C |
| Carrier Gas: | He, 39.0 psi |
| Carrier Gas Flow : | 0.370 mL/min |
| Constant Flow: | On |
| Average Linear Velocity: | 35 cm/sec at 70 °C |
| Injection Mode: | Split |
| Purge On Time: | 1.0 |
| Purge On (Split) Vent Flow: | 10 mL/min |
| Injection Volume: | 0.5 µL |
| Injection Temperature: | 240 °C |
| Autosampler: | Yes |
| Liner Type : | 2.3 mm ID FocusLiner™ |
| Liner Part Number: | 092005 |



GC Columns and Applications

NORMAL

Chromatogram showing separation of Organophosphorous Pesticides using a conventional 30 meter x 0.25 mm ID BPX50 column with a 0.25 micron film.




Components

- | | |
|--------------------------------------|----------------------|
| 1. 4-Chloro-3-nitrobenzo-trifluoride | 8. Naled |
| 2. Dichlorvos | 9. Phorate |
| 3. 1-Bromo-2-nitrobenzene | 10. Demeton |
| 4. α-Mevinphos | 11. Diazinon |
| 5. Tributylphosphate (IS) | 12. Disulfoton |
| 6. Ethoprop | 13. Methyl Parathion |
| 7. Sulfotepp | 14. Ronnel |
| | 15. Chlorpyrifos |
| | 16. Fenthion |

FAST

Chromatogram showing separation of Organophosphorous Pesticides using a FAST BPX50 column.



- | |
|-----------------------------|
| 17. Trichlorinate |
| 18. Tetrachlorvinphos |
| 19. Tokuthion |
| 20. Impurity |
| 21. Fensulfothion |
| 22. Impurity |
| 23. Triphenylphosphate (IS) |
| 24. Guthion |
| 25. Coumaphos |





GC Columns and Applications

ENV 04 | Analysis of Herbicides on BPX35

Column Part No.: 054711

Phase: BPX35, 0.25 µm film

Column: 25 m x 0.22 mm ID

Initial Temp.: 80 °C

Rate: 10 °C/min


Final Temp.: 300 °C 5 min

Detector: FID, 380 °C

Injection Mode: Split (20:1)

Carrier Gas: He, 100 kpa

Note: BPX35 provides quick analysis of all 3 Triazine compounds



Components

1. Eptam®
2. Sutan®
3. Vernam®
4. Tillam®
5. Ordram®
6. Treflan®
7. Balan®
8. Ro-Neet®
9. Propachlor
10. Tolban®
11. Propazine
12. Atrazine
13. Simazine
14. Terbacil
15. Sencor®
16. Dual®
17. Paarlan®
18. Prowl®
19. Bromacil
20. Oxadiazon
21. GOAL®
22. Hexazinone

ENV 48 | Analysis of Herbicides on BPX5

Column Part No.: 054101

Phase: BPX5, 0.25 µm

Column: 30 m x 0.25 mm ID

Initial Temp.: 90 °C, 1 min

Rate 1: 30 °C/min

Temp.: 180 °C


Rate 2: 5 °C/min

Final Temp.: 260 °C, 10 min

Detector: NPD


Injection Mode: Varian SPI

Carrier Gas: He, 10 psi



Components

1. Metamidofos
2. Acephate
3. Diphenylamine
4. Monocrofos
5. Sulfotep
6. Tolclofos-methyl
7. Fenitrothion
8. Triadimefon
9. Trichloronate
10. Triadimenol
11. Bromophos-ethyl
12. Bupirimate
13. Carbophenothion
14. Dialifos






GC Columns and Applications

ARO 14 | Analysis of chlorinated and nitroaromatic compounds on SolGel-1ms™

| | | | |
|-------------------------|----------------------------|--------------------------|----------------------------|
| Column Part No.: | 054462 | Constant Flow: | On |
| Phase: | SolGel-1ms™ 0.25 µm film | Average Linear Velocity: | 35 cm/sec, 40 °C |
| Sample: | 200 ppm in dichloromethane | Injection Mode: | Split |
| Column: | 30 m x 0.25 mm ID | Split Ratio: | 100 : 1 |
| Initial Temp: | 40 °C, 1 min. | Injection Volume: | 0.5 µL |
| Rate 1: | 10 °C/min to 300 °C | Injection Tem: | 250 °C |
| Final Temp: | 300 °C, 2 min. | Liner Type: | 4 mm ID Single Taper Liner |
| Detector Type: | MSD | Liner Part No.: | 092017 |
| Carrier Gas: | He, 25.7 psi | Full Scan / SIM: | Full scan 45-450 |
| Carrier Gas Flow: | 1.8 mL/min. | | |




Components

1. Phenol
2. o-Chlorophenol
3. p-Dichlorobenzene
4. Nitrobenzene
5. o-Nitrophenol
6. 2,4-Xylenol
7. 2,4-Dichlorophenol
8. 4-Chloro-3-methylphenol
9. 2,4,6-Trichlorophenol
10. 2,6-Dinitrotoluene
11. 2,4-Dinitrophenol
12. 2,4-Dinitrotoluene
13. 4-Nitrophenol
14. 4,6-Dinitro-o-cresol
15. 4-Chlorophenyl phenyl ether
16. 4-Bromophenyl phenyl ether
17. Hexachlorobenzene
18. Pentachlorophenol

ALC 06 | US EPA 625 Phenols Mix on BPX50

| | | | |
|-------------------------|-------------------|--------------------|----------------|
| Column Part No.: | 054751 | Initial Oven Temp: | 50 °C, 1 min |
| Phase: | BPX50, 0.25 µm | Rate 1: | 8 °C/min |
| Column: | 30 m x 0.25 mm ID | Final Temp: | 300 °C, 10 min |
| Injector Mode: | Split, 40:1 | Detector: | HP 5973 MSD |



Components

1. 2-Chlorophenol
2. 2-Nitrophenol
3. 2, 4-Dimethylphenol
4. 2, 4-Dichlorophenol
5. 4-Chloro-3-methylphenol
6. 2, 4, 6-Trichlorophenol
7. 2, 4- Dinitrophenol
8. 4-Nitrophenol
9. 2-Methyl-4, 6-dinitrophenol
10. Pentachlorophenol






GC Columns and Applications



ENV 65 | Analysis of Synthetic Pyrethroids on BPX50

| | | | |
|-------------------------|----------------------|-----------------------------|----------------------------|
| Column Part No.: | 054751 | | |
| Phase: | BPX50, 0.25 µm film | Constant Flow: | On |
| Column: | 30 m x 0.25 mm ID 16 | Average Linear Velocity: | 36 cm/sec at 50 °C |
| Pyrethroids: | 10 ppm in methanol | Injection Mode: | Splitless |
| Initial Temp.: | 50 °C, 1 min. | Purge on Time: | 0.5 min |
| Rate 1: | 30 °C/min to 200 °C | Purge on (Split) Vent Flow: | 60 mL/min |
| Rate 2: | 4 °C/min to 300 °C | Injection Volume: | 1 µL |
| Final Temp.: | 300 °C, 5 min | Injection Temperature: | 250 °C |
| Detector Type: | MSD | Autosampler: | No |
| Carrier Gas: | He, 6.8 psi | Liner Type: | 4 mm ID Double Taper Liner |
| Carrier Gas Flow: | 1.0 mL/min | Liner Part Number: | 092018 |



ARO 13 | Analysis of BTEX on SolGel-WAX™

| | | | |
|-------------------------|--------------------------|--------------------------|---------------------------|
| Column Part No.: | 054796 | | |
| Phase: | SolGel-WAX™ 0.25 µm film | Constant Flow: | On |
| BTEX: | 300 ppm in methanol | Average Linear Velocity: | 35 cm/sec, 60 °C |
| Column: | 30 m x 0.25 mm ID | Injection Mode: | Split |
| Initial Temp: | 60 °C, 10 min | Split Ratio: | 100:1 |
| Detector Type: | FID | Injection Volume: | 0.2 µL |
| Carrier Gas: | He, 17.3 psi | Injection Temp: | 250 °C |
| Carrier Gas Flow: | 1.5 mL/min | Liner Type: | 4 mm ID Double Taper Line |
| | | Liner Part Number: | 092018 |




ENV 84 | Analysis of US EPA 8270 Mix on BPX5

| | | | |
|-------------------------|--|-----------------------------|--------------------------|
| Column Part No.: | 054101 | | |
| Phase: | BPX5, 0.25 µm film | Carrier Gas Flow: | 1.1 mL/min. |
| Column: | 30 m x 0.25 mm | Constant Flow: | On |
| ID Sample: | 5 ppm solution | Injection Mode: | Splitless |
| Initial Temp.: | 40 °C, 3 min | Purge on Time: | 0.5 min |
| Rate 1: | 8 °C/min to 300 °C | Purge on (Split) Vent Flow: | 40 mL/min |
| Final Temp.: | 300 °C, 9 min. | Injection Volume: | 1 µL |
| Detector Type: | Mass Spectrometer | Injection Temperature: | 250 °C |
| Carrier Gas: | He | Autosampler: | No |
| Inlet Pressure: | 16 psi for 30 sec then drops to 10 psi | Liner Type: | 4 mm ID Single Gooseneck |
| Pressure rate1: | 10 psi to 28 psi at 0.5 psi/min | Liner Part Number: | 092017 |
| Final Pressure: | 28 psi until end of run | Full Scan / SIM: | Full scan 41-450 |



GC Columns and Applications



Components





- | | | |
|-----------------------------------|---------------------------------|----------------------------------|
| 1. Pyridine | 28. Naphthalene-d8 | 57. 2,4,6-Tribromophenol |
| 2. n-Nitrosodimethylamine | 29. Naphthalene | 58. 4-Bromophenyl phenyl ether |
| 3. 2-Fluorophenol | 30. Hexachlorobutadiene | 59. Hexachlorobenzene |
| 4. Phenol-d5 | 31. 4-Chloroaniline | 60. Pentachlorophenol |
| 5. Phenol | 32. 4-Chloro-3-methylphenol | 61. Phenanthrene-d10 |
| 6. Aniline | 33. 2-Methylnaphthalene | 62. Phenanthrene |
| 7. 2-Chlorophenol | 34. Hexachlorocyclopentadiene | 63. Anthracene |
| 8. bis- (2-chloroethyl) ether | 35. 2,4,6-Trichlorophenol | 64. Carbazole |
| 9. 1,3-Dichlorobenzene | 36. 2,4,5-Trichlorophenol | 65. Di-n-butyl phthalate |
| 10. 1,4-Dichlorobenzene-d4 | 37. 2-Fluorobiphenyl | 66. Fluoranthene |
| 11. 1,4-Dichlorobenzene | 38. 2-Chloronaphthalene | 67. Benzidine |
| 12. 1,2-Dichlorobenzene | 39. 2-Nitroaniline | 68. Pyrene |
| 13. Benzyl alcohol | 40. Dimethyl phthalate | 69. p-Terphenyl-d14 |
| 14. 2-Methyl phenol | 41. Acenaphthylene | 70. Butyl benzyl phthalate |
| 15. bis-(2-chloroisopropyl)ether | 42. 2,6-Dinitrotoluene | 71. Benz[a]anthracene |
| 16. n-Nitroso-di-n-propylamine | 43. Acenaphthene-d10 | 72. Chrysene-d12 |
| 17. Hexachloroethane | 44. Acenaphthene | 73. Chrysene |
| 18. 4-Methylphenol | 45. 3-Nitroaniline | 74. 3,3-Dichlorobenzidine |
| 19. Nitrobenzene-d5 | 46. 2,4-Dinitrophenol | 75. bis (2-Ethylhexyl) phthalate |
| 20. Nitrobenzene | 47. Dibenzofuran | 76. Di-n-octyl phthalate |
| 21. Isophorone | 48. 4-Nitrophenol | 77. Benzo (b) fluoranthene |
| 22. 2-Nitrophenol | 49. 2,4-Dinitrotoluene | 78. Benzo (k) fluoranthene |
| 23. 2,4-Xylenol | 50. Diethylphthalate | 79. Benzo (a) pyrene |
| 24. bis- (2-Chloroethoxy) methane | 51. Fluorene | 80. Perylene-d12 |
| 25. Benzoic acid | 52. 4-Chlorophenyl phenyl ether | 81. Indeno (1,2,3-cd) perylene |
| 26. 2,4-Dichlorophenol | 53. 2-Methyl-4,6-dinitrophenol | 82. Dibenz (a,h) anthracene |
| 27. 1,2,4-Trichlorobenzene | 54. 4-Nitroaniline | 83. Benzo[g,h,i]perylene |
| | 55. n-Nitrosodiphenylamine | |
| | 56. Azobenzene | |





GC Columns and Applications

ENV 84 continued




SGE would like to thank Mark Ferry from ECS/MDL USA for supplying all of the chromatograms for this application note.

FOO 03 | Analysis of Scotch Whisky on BP20

| | |
|-------------------------|-----------------------------|
| Column Part No.: | 054447 |
| Phase: | BP20, 1.0 µm film |
| Column: | 12 m x 0.53 mm ID |
| Initial Temp: | 55 °C, 3 min |
| Rate: | 10 °C/min |
| Final Temp: | 120 °C, 0 min |
| Detector: | FID |
| Sensitivity: | 128 x 10 ⁻¹² AFS |
| Injection Mode: | Split |

Components



1. Acetaldehyde
2. Ethyl Acetate
3. Methanol
4. Ethanol
5. Propan-1-ol
6. 2-Methylpropan-1-ol
7. 2-Methylbutan-1-ol + 3-Methylbutan-1-ol



GC Columns and Applications

AN-0022-C | FAME Analysis with BPX90 – A Highly Polar Column

| | | | |
|-------------------------|---------------------------------------|-------------------|---------------------------|
| Column Part No.: | 054570 | Constant Flow: | ON |
| Phase: | 90% Cyanopropyl Polysilphenylsiloxane | Pressure: | 4.02 psi |
| Column Dimensions: | 15 m x 0.25 mm x 0.25 µm | Column Flow Rate: | 1.3 ml/min |
| Injector Temperature: | 250 °C | Linear Velocity: | 59 cm/sec |
| Injection Volume: | 1.0 µL | Initial Temp.: | 70 °C hold for 1 minute |
| Injector Type: | Split | Rate: | 20 °C/min to 150 °C |
| Split Ratio: | 100:1 | Rate: | 10 °C/min |
| Liner Type: | FocusLiner™ | Final Temp.: | 250 °C hold for 5 minutes |
| Carrier Gas: | Helium | Detector Type: | MSD |



Supelco 37 FAME standard analyzed with (a) BPX70 and (b) BPX90

SGE would like to thank J. Harynyk, P.J. Marriott and P. Wynne, *Chromatographia*, 2006; 63 (Supplement 13): S61-S66.






GC Columns and Applications

AN-0011-C | Analysis of Omega-3 Fatty Acids using a Highly Selective GC Capillary Column

| | |
|-------------------------|-----------------------------|
| Column Part No.: | 054606 |
| Phase: | BPX70, 0.25 µm film |
| Sample: | 10 ppm in methanol |
| Column: | 25 m x 0.32 mm ID |
| Initial Temp: | 80 °C, 2 min |
| Rate 1: | 50 °C/min to 130 °C, 10 min |
| Rate 2: | 2 °C/min to 172 °C |
| Final Temp: | 172 °C, 6 min |
| Detector Type: | FID |
| Detector Temp: | 300 °C |
| Carrier Gas: | He, 10 psi |

| | |
|--------------------------|---------------------|
| Carrier Gas Flow: | 2.2 mL/min |
| Constant Flow: | On |
| Average Linear Velocity: | 39 cm/sec at 80 °C |
| Injection Mode: | Split |
| Split Ratio: | 58:1 |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Autosampler: | No |
| Liner Type: | 4 mm ID FocusLiner™ |
| Liner Part No.: | 092002 |



Components:

1. C6
2. C8
3. C10
4. C11
5. C12
6. C13
7. C14
8. C14:1 cis 9
9. C15
10. C16
11. C16:1 cis 9

12. C17
13. C17:1 cis 10
14. C18
15. C18:1 trans 9
16. C18:1 cis 9
17. C18:1 cis 12
18. C18:2 trans 9, 12
19. C18:2 cis 9, 12
20. C18:3 cis 6, 9, 12
21. C18:3 cis 9, 12, 15
22. C18:4 Cis 6, 9, 12, 15
23. C20


24. C21:1 cis 11
25. C20:2 cis 11,14
26. C20:3 cis 8, 11,14
27. C20:4 cis 5, 8, 11, 14
28. C20:3 cis 11, 14, 17
29. C22
30. C20:5 cis 5, 8, 11, 14, 17
31. C22:1 cis 13
32. C22:4 cis 7, 10, 13, 16
33. C24
34. C22:6 cis 4,7,10,13,16,19

Notes: The chromatogram shows the excellent separation of a complex mixture of FAME compounds. Note the excellent peak shape and separation of the Omega-1,2 and 3 fatty acid isomers both structural and cis and trans.


SGE would like to thank Masterfoods UK for supplying the sample and chromatographic conditions for this chromatogram.

FOO 16 | Analysis of Triglyceride Standards on HT5

| | |
|-------------------------|------------------------------------|
| Column Part No.: | 054661 |
| Phase: | HT5, 0.1 µm |
| Column: | 6 m x 0.53 mm I.D. (Aluminum Clad) |
| Initial Temp.: | 60 °C, 0 min |
| Program Rate: | 10 °C/min |
| Final Temp.: | 370 °C, 5 min |
| Carrier Gas: | H ₂ , 2 psi |
| Detector: | F.I.D. |
| Sensitivity: | 32 x 10 ⁻¹² AFS |
| Injection Mode: | On-column |



Notes: For the analysis of triglycerides, on-column injection is recommended. Temperatures above 380 °C are not recommended as triglycerides can degrade.




FLA 05 | Analysis of Menthol Oil on CYDEX-B

| | | | |
|-------------------------|-----------------------|-----------------|----------------|
| Column Part No.: | 054901 | Final Temp.: | 130 °C |
| Phase: | Cydex-B, 0.25 µm film | Carrier Gas: | H ₂ |
| Column: | 50 m x 0.22 mm I.D. | Detector: | F.I.D. |
| Initial Temp.: | 100 °C, 5 min | Sensitivity: | 32 x 10-12 AFS |
| Rate: | 2 °C/min | Injection Mode: | Split |



GC Columns and Applications




Components

1. (+) Neomenthol
2. (-) Neomenthol
3. (+) Menthol
4. (-) Menthol
5. (+) a-Terpineol
6. (-) a-Terpineol

Notes: Cydex - B column enables the separation of three different enantiomer pairs in Menthol Oil.

FLA 19 | Analysis of Eucalyptus Oil on SolGel-WAX™

| | | | |
|-------------------------|---------------------------|--------------------------|----------------------------|
| Column Part No.: | 054796 | Constant Flow: | On |
| Phase: | SolGel-WAX™, 0.25 µm film | Average Linear Velocity: | 35 cm/sec at 40 °C |
| Sample: | Neat | Injection Mode: | Split |
| Column: | 30 m x 0.25 mm ID | Split Ratio: | 100:1 |
| Initial Temp.: | 40 °C, 1 min. | Injection Volume: | 0.2 µL |
| Rate 1: | 8 °C/min to 220 °C, | Injection Temp.: | 250 °C |
| Final Temp: | 220 °C, 5 min. | Liner Type: | 4 mm ID Single Taper Liner |
| Detector Type: | Mass Spectrometer | Liner Part Number: | 092017 |
| Carrier Gas: | He, 25.7 psi | Full Scan / SIM: | Full scan 45-450 |
| Carrier Gas Flow: | 1.8 mL/min. | | |



Components


1. α-Pinene
2. β-Pinene
3. Sabinene
4. Limonene
5. 1,8-Cineole
6. p-Cymene
7. Linalool L
8. Terpinen-4-ol
9. α-Terpineol
10. α-Terpinenyl acetate
11. d-Carvone





GC Columns and Applications

FOO 02 | Analysis of Food Additives Antimicrobials on BP5



| | |
|-------------------------|-------------------|
| Column Part No.: | 054186 |
| Phase: | BP5, 0.5 µm film |
| Column: | 25 m x 0.32 mm ID |
| Initial Temp: | 160 °C, 0 min |
| Rate: | 15 °C/min |
| Final Temp: | 280 °C, 0 min |
| Detector: | FID |
| Sensitivity: | 256 x 10-12 AFS |
| Injection Mode: | Split |

Components



1. Methyl Paraben
2. Ethyl Paraben
3. Propyl Paraben
4. Butyl Paraben
5. Heptyl Paraben

FOO 04 | Analysis Of 13 Sugar Component Alditol Acetate Mixture on BPX70

| | |
|-------------------------|--|
| Column Part No.: | 054622 |
| Phase: | BPX70, 0.25 µm film |
| Column: | 30 m x 0.25 mm I.D. |
| Initial Temp.: | 190 °C, 1 min. |
| Program Rate: | 3 °C/min. |
| Final Temp: | 260 °C, 10min. |
| Carrier Gas: | He, 50 kPa |
| Detector: | MS (Electron Impact Ionisation, Ionisation Potential 70 eV) Scan 100 m/z to 350 m/z in 0.3 s. |
| Injection Mode: | Split 50:1 |

Components

1. Erythritol
2. 2-Deoxy-ribitol
3. Rhamnitol
4. Fucitol
5. Ribitol
6. Arabinitol
7. Xylitol
8. 2-Deoxy-glucitol
9. Allitol
10. Mannitol
11. Galacitol
12. Glucitol
13. Myo-inositol




FOO 21 | Plant Sterols

| | |
|-------------------------|-------------------|
| Column Part No.: | 054148 |
| Phase: | BPX5, 1.0 µm |
| Column: | 30 m x 0.53 mm ID |
| Initial Temp.: | 320 °C |
| Detector: | FID, 360 °C |
| Injector Mode: | split 100:1 |
| Carrier Gas: | He, 3 psi |
| Injection Volume: | 1 µL |

Components

- | | |
|-------------------|---------------------|
| 1. Coprosterol | 7. Campesterol |
| 2. Cholesterol | 8. Stigmasterol |
| 3. Cholestanol | 9. Unknown |
| 4. Desmosterol | 10. beta-Sitosterol |
| 5. Brassicasterol | 11. Lanosterol |
| 6. Ergosterol | |




GC Columns and Applications

FOO 24 | Analysis of Wine on BP20

| | |
|-------------------------|------------------------|
| Column Part No.: | 054442 |
| Phase: | BP20, 1.0 µm |
| Column: | 25 m x 0.32 mm ID |
| Initial Temp.: | 40 °C, 2 min |
| Rate 1: | 5 °C/min |
| Temp 2: | 50 °C |
| Rate 2: | 15 °C/min |
| Final Temp.: | 190 °C |
| Carrier Gas: | H ₂ , 6 psi |
| Injection Mode: | 2 µL |

Components

- | | |
|------------------|--------------------|
| 1. Acetaldehyde | 6. Propanol |
| 2. Ethyl Acetate | 7. Isobutanol |
| 3. Methanol | 8. Isoamyl Alcohol |
| 4. Isopropanol | 9. Acetic Acid |
| 5. Ethanol | 10. Unknown |






GC Columns and Applications

FLA 14 | Analysis of Lavender Oil on BPX5



| | |
|-------------------------|---------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.2 5 µm film |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp.: | 40 °C, 1 min |
| Rate 1: | 5 °C/min to 260 °C |
| Final Temp.: | 260 °C |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 7.0 psi |
| Carrier Gas Flow: | 1.0 mL/min |
| Constant Flow: | On |

| | |
|--------------------------|----------------------------|
| Average Linear Velocity: | 36 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 200:1 |
| Purge on (Split) | |
| Vent Flow: | 200 mL/min |
| Injection Volume: | 0.2 µL |
| Injection Temp.: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |



FLA 03 | Analysis of Orris Concentrate on BP20

| | |
|-------------------------|---------------------|
| Column Part No.: | 054436 |
| Column: | BP20, 0.5 µm |
| Phase: | 25 m x 0.32 mm I.D. |
| Initial Temp.: | 70 °C, 1 min |
| Rate: | 10 °C/min |
| Final Temp.: | 250 °C, 10 min |
| Carrier Gas: | Helium |
| Carrier Pressure: | 10 psi |
| Injection Mode: | Split 50:1 |




FLA 21 | Analysis of Spearmint Oil on SolGel-WAX™



GC Columns and Applications

| | | | |
|--------------------------|---------------------------|---------------------------------|----------------------------|
| Column Part No.: | 054796 | Constant Flow: | On |
| Phase: | SolGel-WAX™, 0.25 µm film | Average Linear Velocity: | 35 cm/sec at 40 °C |
| Sample: | Neat | Injection Mode: | Split |
| Column: | 30 m x 0.25 mm ID | Split Ratio: | 100:1 |
| Initial Temp.: | 40 °C, 1 min. | Injection Volume: | 0.2 µL |
| Rate 1: | 8 °C/min to 220 °C | Injection Temp.: | 250 °C |
| Final Temp: | 220 °C, 5 min. | Liner Type: | 4 mm ID Single Taper Liner |
| Detector Type: | Mass Spectrometer | Liner Part Number: | 092017 |
| Carrier Gas: | He, 25.7 psi | Full Scan / SIM: | Full scan 45-450 |
| Carrier Gas Flow: | 1.8 mL/min. | | |




Components

- | | |
|------------------------|---------------------------------|
| 1. α -Pinene | 11. <i>trans</i> Caryophyllene |
| 2. β -Pinene | 12. <i>cis</i> dihydrocarvone |
| 3. Sabinene | 13. <i>trans</i> dihydrocarvone |
| 4. Myrcene | 14. Menthol |
| 5. Limonene | 15. Dihydrocarvyl acetate |
| 6. 1,8-Cineole | 16. L-Carvone |
| 7. 3-Octanol | 17. <i>trans</i> Carveol |
| 8. Menthone | 18. <i>cis</i> Carveol |
| 9. β -bourbonene | 19. Caryophyllene oxide |
| 10. Linalool L | |

FLA 18 | Analysis of Ylang Ylang Oil on SolGel-WAX™

| | | | |
|--------------------------|-----------------------|---------------------------------|----------------------------|
| Column Part No.: | 054796 | Constant Flow: | On |
| Phase: | SolGel-WAX™, 0.25 µm | Average Linear Velocity: | 35 cm/sec at 40 °C |
| Sample: | Ylang Ylang oil neat. | Injection Mode: | Split |
| Column: | 30 m x 0.25 mm ID | Split Ratio: | 120:1 |
| Initial Temp.: | 40 °C, 2 min. | Injection Volume: | 0.1 µL |
| Rate 1: | 3 °C/min to 250 °C | Injection Temp.: | 250 °C |
| Final Temp: | 250 °C, 10 min. | Autosampler: | No |
| Detector Type: | Mass Spectrometer | Liner Type: | 4 mm ID Double Taper Liner |
| Carrier Gas: | He, 25.7 psi | Liner Part Number: | 092018 |
| Carrier Gas Flow: | 1.8 mL/min. | Full Scan / SIM: | Full scan 45-450 |



Components

1. p-Methyl anisole
2. α -Copaene
3. Linalool L
4. β -Caryophyllene
5. Methyl benzoate
6. α -Humulene
7. α -Amorphene
8. Germacrene
9. Benzyl acetate
10. δ -Cadinene
11. α -Farnesene
12. Geranyl acetate
13. *trans*-Cinamyl acetate
14. Farnesyl acetate
15. Farnesol
16. Benzyl benzoate
17. Benzyl salicylate






GC Columns and Applications

FLA 16 | Analysis of Pine Oil on BPX5

| | |
|-------------------------|--------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp.: | 40 °C, 1 min. |
| Rate 1: | 5 °C/min to 260 °C |
| Final Temp: | 260 °C |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 7.0 psi |
| Carrier Gas Flow: | 1.0 mL/min. |
| Constant Flow: | On |

| | |
|--------------------------|----------------------------|
| Average Linear Velocity: | 36 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 200:1 |
| Purge on (Split) | |
| Vent Flow: | 200 mL/min. |
| Injection Volume: | 0.2 µL |
| Injection Temp.: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |




Components

1. α-Pinene
2. Camphene
3. β-Pinene
4. δ-3-Carene
5. Limonene
6. Endobornyl acetate

FLA 15 | Analysis of Tea Tree Oil on BPX5


| | |
|-------------------------|--------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp.: | 40 °C, 1 min. |
| Rate 1: | 5 °C/min to 200 °C |
| Final Temp: | 200 °C |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 7.0 psi |
| Carrier Gas Flow: | 1.0 mL/min. |
| Constant Flow: | On |

| | |
|--------------------------|----------------------------|
| Average Linear Velocity: | 36 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 200:1 |
| Purge on (Split) | |
| Vent Flow: | 200 mL/min. |
| Injection Volume: | 0.2 µL |
| Injection Temp.: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |



Components

- | | |
|-------------------|-----------------------------|
| 1. Thujene | 12. Terpinolene |
| 2. α-Pinene | 13. Terpinen-4-ol |
| 3. Sabinene | 14. α-Terpineol |
| 4. 3-Octanol | 15. α-Gurjunene |
| 5. β-Pinene | 16. (trans)-β-Caryophyllene |
| 6. α-Phellandrene | 17. Aromadendrene |
| 7. α-Terpinene | 18. Alloaromadendrene |
| 8. p-Cymene | 19. Ledene |
| 9. Limonene | 20. Germacrene B |
| 10. 1,8-Cineole | 21. δ-Cadinene |
| 11. γ-Terpinene | 22. 1s, cis-Calamenene |
| | 23. Globulol |



FLA 12 | Analysis of Nutmeg Oil on BPX5

| | |
|-------------------------|---------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp.: | 40 °C, 1 min. |
| Rate: | 5 °C/min to 260 °C, |
| Final Temp: | 260 °C |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 7.0 psi |
| Carrier Gas Flow: | 1.0 mL/min. |


| | |
|-----------------------------|----------------------------|
| Constant Flow: | On |
| Average Linear Velocity: | 36 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 200:1 |
| Purge on (Split) Vent Flow: | 200 mL/min. |
| Injection Volume: | 0.2 µL |
| Injection Temp.: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |



GC Columns and Applications

Components

- | | |
|-------------------|----------------------------|
| 1. α-Thujene | 11. Limonene |
| 2. α-Pinene | 12. β-Phellandrene |
| 3. Camphene | 13. γ-Terpinene |
| 4. Sabinene | 14. trans-Sabinene hydrate |
| 5. β-Pinene | 15. α-Terpinolene |
| 6. Myrcene | 16. Terpinen-4-ol |
| 7. α-Phellandrene | 17. 1-α-Terpineol |
| 8. δ-3-Carene | 18. Safrole |
| 9. α-Terpinene | 19. α-Copaene |
| 10. p-Cymene | 20. Myristicin |
| | 21. Elemicin |






GC Columns and Applications

FLA 13 | Analysis of Rosemary Oil on BPX5


| | |
|-------------------------|---------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp.: | 40 °C, 1 min. |
| Rate 1: | 5 °C/min to 260 °C, |
| Final Temp: | 260 °C |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 7.0 psi |
| Carrier Gas Flow: | 1.0 mL/min. |
| Constant Flow: | On |

| | |
|--------------------------|----------------------------|
| Average Linear Velocity: | 36 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 200:1 |
| Purge on (Split) | |
| Vent Flow: | 200 mL/min. |
| Injection Volume: | 0.2 µL |
| Injection Temp.: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |



Components

- | | |
|-------------------------|-------------------------|
| 1. α-Pinene | 12. Camphor |
| 2. Camphene | 13. Borneol L |
| 3. Sabinene | 14. Terpinen-4-ol |
| 4. β-Pinene | 15. α-Terpineol |
| 5. Myrcene | 16. Verbenone |
| 6. p-Cymene | 17. trans-(+)-Carveol |
| 7. Limonene | 18. Bornyl acetate |
| 8. 1,8-Cineole | 19. α-Copaene |
| 9. cis-Sabinene hydrate | 20. β-Caryophyllene |
| 10. Linalool L | 21. α-Humulene |
| 11. Isopinocarveol | 22. α-Amorphene |
| | 23. Caryophyllene oxide |




FLA 22 | Analysis of Tasmanian Lavender Oil on SolGel-WAX™



GC Columns and Applications

| | | | |
|-------------------------|---------------------------|--------------------------|----------------------------|
| Column Part No.: | 054796 | | |
| Phase: | SolGel-WAX™, 0.25 µm film | Constant Flow: | On |
| Sample: | Neat | Average Linear Velocity: | 35 cm/sec at 40 °C |
| Column: | 30 m x 0.25 mm ID | Injection Mode: | Split |
| Initial Temp.: | 40 °C, 1 min. | Split Ratio: | 100:1 |
| Rate 1: | 8 °C/min to 220 °C, | Injection Volume: | 0.2 µL |
| Final Temp: | 220 °C, 5 min. | Injection Temp.: | 250 °C |
| Detector Type: | Mass Spectrometer | Liner Type: | 4 mm ID Single Taper Liner |
| Carrier Gas: | He, 25.7 psi | Liner Part Number: | 092017 |
| Carrier Gas Flow: | 1.8 mL/min. | Full Scan / SIM: | Full scan 45-450 |



Components

1. 3-Octanone
2. Octenyl acetate
3. Octanol
4. Cis Linalool oxide
5. Trans Linalool Oxide
6. Linalool L
7. Linalyl acetate
8. Terpinen-4-ol
9. Lavandulyl acetate
10. Borneol L
11. Caryophyllene oxide





GC Columns and Applications


ENV 54 | Total Recoverable Petroleum Hydrocarbons (TRPH) Analysis on Standard and Fast BPX5

| | |
|-----------------------------|---------------------------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| TRPH (C8-C40): | 5 ng/ µL in dichloromethane |
| Initial Temp: | 40 °C , 2 min |
| Rate 1: | 30 °C/min to 330 °C |
| Rate 2: | N/A |
| Final Temp.: | 330 °C, 9 min |
| Detector Type: | FID, 350 °C |
| Carrier Gas: | He, 14.1 psi |
| Carrier Gas Flow : | 1.29 mL/min |
| Constant Flow: | On |
| Average Linear Velocity: | 40 cm/sec at 40 °C |
| Injection Mode: | Split, 120:1 |
| Purge On Time: | N/A |
| Purge On (Split) Vent Flow: | 160 mL/min |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Autosampler: | Yes |
| Liner Type : | 4 mm ID FocusLiner™ with single taper |
| Liner Part Number: | 092003 |

| | |
|-----------------------------|-----------------------------|
| Column Part Number: | 054099 |
| Phase: | BPX5, 0.10 µm film |
| Column: | 10 m x 0.10 mm ID |
| TRPH (C8-C40) Standard: | 5 ng/ µL in dichloromethane |
| Initial Temp.: | 40 °C , 1 min |
| Rate 1: | 30 °C/min to 330 °C |
| Rate 2: | N/A |
| Final Temp: | 330 °C, 0 min |
| Detector Type: | FID, 350 °C |
| Carrier Gas: | He, 28 psi |
| Carrier Gas Flow : | 0.52 mL/min |
| Constant Flow: | On |
| Average Linear Velocity: | 55 cm/sec at 40 °C |
| Injection Mode: | Split, 120:1 |
| Purge On Time: | N/A |
| Purge On (Split) Vent Flow: | 62 mL/min |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Autosampler: | Yes |
| Liner Type : | 2.3 mm ID FocusLiner™ |
| Liner Part Number: | 092005 |



NORMAL

Chromatogram showing separation of Total Recoverable Petroleum Hydrocarbons using a conventional 30 meter x 0.25 mm ID BPX5 column with a 0.25 micron film.



FAST

Chromatogram showing separation of Total Recoverable Petroleum Hydrocarbon using a FAST BPX5 column.




PET 01 | Analysis of Gasoline Range Hydrocarbons on BP1-PONA



GC Columns and Applications

| | | | |
|-------------------------|-------------------|-----------------|-------------------------|
| Column Part No.: | 054950 | Temp. 3: | 120 °C |
| Phase: | BP1 PONA | Rate 3: | 10 °C/min |
| Column: | 50 m x 0.15 mm ID | Final Temp.: | 190 °C |
| Initial Temp.: | 30 °C, 5 min hold | Detector: | FID |
| Rate 1: | 2 °C/min | Sensitivity: | 32 x 10-12 AFS |
| Temp. 2: | 80 °C | Injection Mode: | Split |
| Rate 2: | 50 °C/min | Carrier Gas: | H ₂ , 40 psi |



Components

| TIME | COMPOUND |
|-------|---|
| 4.85 | Cyclopentane |
| 5.00 | 2,3-Dimethylbutane |
| 5.25 | 2-Methylpentane |
| 5.74 | 3-Methylpentane |
| 6.45 | n-Hexane |
| 7.46 | 2,2-Dimethylpentane |
| 7.60 | Methylcyclopentane |
| 7.91 | 2,4-Dimethylpentane |
| 8.18 | 2,2,3-Trimethylbutane |
| 8.99 | Benzene |
| 9.35 | 3,3-Dimethylpentane |
| 9.55 | Cyclohexane |
| 10.23 | 2-Methylhexane |
| 10.32 | 2,3-Dimethylpentane |
| 10.47 | 1,1-Dimethylcyclohexane |
| 10.83 | 3-Methylhexane |
| 11.23 | 1-trans-3-Dimethylcyclopentane |
| 11.43 | 1-cis-3-Dimethylcyclopentane |
| 11.55 | 3-Ethylpentane |
| 11.63 | 1-trans-2-Dimethylcyclopentane |
| 11.78 | 2,2,4-Trimethylpentane |
| 12.73 | n-Heptane |
| 14.23 | Methylcyclohexane |
| 14.53 | 2,2-Dimethylhexane |
| 15.27 | Ethylcyclopentane |
| 15.49 | 2,5-Dimethylhexane |
| 15.65 | 2,4-Dimethylhexane |
| 16.09 | 1-trans-2-cis-4-Trimethylcyclopentane |
| 16.24 | 2,3,4-Trimethylpentane |
| 16.78 | 1-trans-2-cis-3-Trimethylcyclopentane |
| 17.05 | 2,3,3-Trimethylpentane |
| 17.39 | Toluene |
| 18.27 | 2,3-Dimethylhexane |
| 18.43 | 2-Methyl-3-ethylpentane |
| 18.84 | 2-Methylheptane |
| 19.69 | 1-Methyl-2-ethylcyclopentane |
| 18.98 | 4-Methylheptane |
| 19.23 | 1-cis-2-cis-4-trans-Trimethylcyclopentane |
| 19.50 | 3-Methylheptane |
| 19.77 | 1-trans-4-Dimethylcyclohexane |
| 20.73 | 1-Methyl-cis-2-ethylcyclopentane |
| 20.86 | 1-Methyl-trans-3-ethylcyclopentane |
| 21.08 | 1-Methyl-cis-3-ethylcyclohexane |
| 21.27 | 1-Ethyl-1-methylcyclopentane |
| 21.53 | 1-trans-2-Dimethylcyclohexane |
| 22.43 | n-Octane |
| 23.05 | iso-Propylcyclopentane |
| 24.14 | 2.2.5-Trimethylhexane |
| 24.19 | 2,2,4-Trimethylhexane |
| 24.53 | 2,4,4-Trimethylhexane |
| 24.79 | 2,3,5-Trimethylhexane |
| 25.16 | 2,4-Dimethylheptane |
| 25.41 | n-Propylcyclopentane |
| 25.73 | 1-cis-2-Dimethylcyclohexane |
| 26.00 | 1,1,3-Trimethylcyclohexane |
| 26.25 | 2,5-Dimethylheptane |
| 26.44 | 3,3-Dimethylheptane |
| 26.58 | 3,5-Dimethylheptane |
| 26.77 | 4,4-Dimethylheptane |
| 26.94 | 2,3,3-Trimethylhexane |
| 27.43 | Ethylbenzene |
| 27.57 | 1-cis-3-cis-5-Trimethylpentane |
| 27.69 | 1,1,4-Trimethylcyclohexane |
| 27.88 | 2,3,4-Trimethylhexane |
| 28.15 | 3,3,4-Trimethylhexane |
| 28.42 | m-Xylene |
| 28.54 | p-Xylene |
| 28.74 | 2,3-Dimethylheptane |
| 28.84 | 1-cis-2-trans-4-trans-Trimethylcyclohexane |
| 28.95 | 1-cis-2-trans-4-cis-Trimethylcyclohexane |
| 29.16 | 3,4-Dimethylheptane |
| 29.31 | 3-Methylethylhexane |
| 29.68 | 4-Methyloctane |
| 29.81 | 2-Methyloctane |
| 30.56 | 3-Methyloctane |
| 30.93 | o-Xylene |
| 31.75 | 1-Methyl-2-propylcyclopentane and 1-Methyl-trans-4-ethylcyclohexane |
| 31.98 | 1-Methyl-cis-4-ethylcyclohexane |
| 32.46 | 3,3-Diethylpentane |
| 32.89 | 2,2,6-Trimethylheptane |
| 33.17 | 1,1,2-Trimethylcyclohexane |
| 33.52 | n-Nonane |
| 34.26 | iso-Propylbenzene |
| 34.48 | tert-Butylcyclopentane |
| 34.68 | tert-Butylbenzene |
| 35.57 | sec-Butylcyclopentane |
| 36.33 | 3-Methylnonane |
| 36.56 | n-Propylbenzene |
| 36.83 | n-Propylcyclohexane |
| 37.12 | m-Ethyltoluene |
| 37.24 | p-Ethyltoluene |
| 37.64 | 1,3,5-Trimethylbenzene |
| 38.20 | 2-Methylnonane |
| 38.36 | o-Ethyltoluene |
| 38.75 | 3,6-Dimethyloctane |
| 38.75 | 1,2,4-Trimethylbenzene |
| 40.32 | n-Decane |
| 40.63 | 1,2,3-Trimethylbenzene |
| 41.57 | 4-Methyldecane |
| 41.94 | sec-Butylbenzene |
| 42.45 | n-Butylbenzene |
| 44.54 | n-Undecane |






GC Columns and Applications

PET 11 | Analysis of Crude Oil and Wax Mixtures on HT5

| | | | |
|-------------------------|-------------------------------------|-----------------|----------------------------|
| Column Part No.: | 054635 | Final Temp.: | 480 °C |
| Phase: | HT5, 0.1 µm | Carrier Gas: | H ₂ , 15 psi |
| Column: | 12 m x 0.22 mm I.D. (Aluminum Clad) | Detector: | F.I.D. |
| Initial Temp.: | 35 °C | Sensitivity: | 32 x 10 ⁻¹² AFS |
| Program Rate: | 10 °C/min. | Injection Mode: | Split |

Notes: HT5 is the best column for the analysis of hydrocarbons C₁₀ - C₇₀.





PET 26 | Standard for D2887 on BPX1

| | | | |
|-------------------------|--------------------|-----------------|---------------|
| Column Part No.: | 054802 | Final Temp.: | 350 °C, 10min |
| Phase: | BPX1, 2.65 µm film | Detector Temp.: | 400 °C |
| Column: | 10 m x 0.53 mm ID | Carrier Gas: | He, 20 mL/min |
| Initial Temp.: | 40 °C | Instrument: | HP 6890 |
| Rate: | 15 °C/min | | |

Separation Systems Injector

| | |
|----------------|----------------|
| Initial Temp.: | 80 °C |
| Rate: | 15 °C/min |
| Final Temp.: | 350 °C, 10 min |




ENV 51 | Total Recoverable Petroleum Hydrocarbons (TRPH) C8-C40 on SolGel-1ms™

| | | | |
|-------------------------|---|------------------|--------------------|
| Column Part No.: | 054795 | | |
| Phase: | SolGel-1, 0.25 µm film 30 m x 0.25 mm ID | Pressure: | 16.6 psi |
| Sample Introduction: | Split / Splitless | Column Flow: | 1.6 mL/min |
| Injector Temp.: | 250 °C | Linear Velocity: | 35 cm/sec at 40 °C |
| Injection Volume: | 0.5 µL | Initial Temp: | 40 °C |
| Autosampler Syringe: | 5 µL Fixed Needle Part No. 001810 | Initial Time: | 2 min |
| Septa: | Auto-Sep T™ Part No. 041882 | Rate 1: | 30 °C/min |
| Injection Type: | Split | Final Temp.1: | 310 °C |
| Purge On Time: | NA | Hold Time: | 0 min |
| Purge On (Spilt) Vent: | 100 mL/min | Rate 2: | 10 °C/min |
| Split Ratio: | 62.5 to 1 | Final Temp. 2: | 340 °C |
| Liner Type: | Double taper Part No. 092018 | Hold Time: | 0 min |
| Carrier Gas: | He | Run Time: | 22.00 min |
| Constant Flow: | On | Detector Type: | FID at 340 °C |


Sample Description: TRPH mix 500 mg/L, 4 ng per component on column.

Components

1. C8
2. C10
3. C12
4. C14
5. C16
6. C18
7. C20
8. C22
9. C24
10. C26
11. C28
12. C30
13. C32
14. C34
15. C36
16. C38
17. C40



PET 27 | Analysis of Polywax 1000 on an Aluminum Clad HT5



| | |
|-------------------------|--------------------|
| Column Part No.: | 054673 |
| Phase: | HT5, 0.075 µm film |
| Column: | 5 m x 0.53 mm ID |
| Initial Temp.: | 40 °C, 1 min |
| Rate: | 10 °C/min |
| Final Temp.: | 450 °C, 10 min |
| Carrier Gas: | He, 20 mL/min |
| Solvent: | CS ₂ |




GC Columns and Applications






GC Columns and Applications



PET 22 | Unleaded Gasoline on BPX5

| | |
|-------------------------|--|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Sample Introduction: | Split / Splitless |
| Injector Temp.: | 240 °C |
| Injection Volume: | 0.1 µL |
| Autosampler Syringe: | 0.5 µL Removable Needle Part No. 000410 |
| Septa: | Auto-Sep T™ Part No. 041882 |
| Injection Type: | Split |
| Purge On Time: | NA |
| Purge On (Split) Vent: | 200 mL/min |
| Split Ratio: | 149 to 1 |
| Liner Type: | FocusLiner™ single taper Part No. 092003 |
| Carrier Gas: | He |

| | |
|------------------|--------------------|
| Constant Flow: | On |
| Pressure: | 13.6 psi |
| Column Flow: | 1.34 mL/min |
| Linear Velocity: | 30 cm/sec at 25 °C |
| Initial Temp.: | 25 °C |
| Initial Time: | 1 min |
| Rate 1: | 30 °C/min |
| Final Temp. 1: | 240 °C |
| Hold Time: | 1 min |
| Run Time: | 9.17 min |
| Final Temp. 2: | 340 °C |
| Hold Time: | 0 min |
| Run Time: | 22.00 min |
| Detector Type: | FID at 280 °C |



Sample Description: Pure Unleaded Gasoline


Components

1. Benzene
2. Toluene
3. Ethylbenzene
4. m, p - Xylene
5. o - Xylene

PET 18 | Analysis of Polywax 655 on Megabore BPX1

| | |
|-------------------------|------------------|
| Column Part No.: | 054800 |
| Phase: | BPX1, 0.1 µm |
| Column: | 5 m x 0.53 mm ID |
| Initial Temp: | 40 °C |
| Rate: | 15 °C |

| | |
|----------------|-----------------|
| Final Temp: | 420 °C, 5 min |
| Detector Temp: | 440 °C |
| Carrier: | He, 10 mL/min |
| Instrument: | HP 6890 |
| Solvent: | CS ₂ |





PET 06 | Analysis of Polywax 655 and Refinery Lubrication Oil on HT5



GC Columns and Applications

| | |
|-------------------------|------------------|
| Column Part No.: | 054661 |
| Phase: | HT5, 0.1 µm |
| Column: | 6 m x 0.53 mm ID |
| Initial Temp.: | 50 °C |
| Rate: | 10 °C/min |
| Final Temp.: | 480 °C, 15 min |

| | |
|-----------------|----------------------------|
| Detector: | FID |
| Sensitivity: | 40 x 10 ⁻¹² AFS |
| Injection Mode: | On-Column |
| Carrier Gas: | Hydrogen, 20 ml/min |
| Solvent: | CS ₂ |



ENV 54 | BPX1 A New Era in Simulated Distillation Technology (SimD)

| | |
|------------------------|------------------|
| Column Part No: | 054800 |
| Phase: | BPX1, 0.1 µm |
| Column: | 5 m x 0.53 mm ID |
| Initial Temp.: | 40 °C |
| Rate: | 15 °C |
| Final Temp.: | 420 °C, 5 min. |

| | |
|----------------|-------------------|
| Detector Temp: | 440 °C |
| Carrier Gas: | Helium, 10 mL/min |
| Instrument: | HP6890 |
| Initial Temp.: | 40 °C |
| Rate: | 15 °C |
| Final Temp.: | 420 °C, 5 min. |

Data supplied by Dr. J. Lubkowitz and the staff at Separation Systems Inc.

Figure. 1. Standard mix for HTSD using BPX1-SimD





Figure. 2. Enlarged section of Figure 1.





GC Columns and Applications

ALC 02 | Analysis of 18 Alcohols on BP20

Column Part No.: 054427

Phase: BP20, 0.25 µm film

Column: 30 m x 0.25 mm ID

Initial Temp: 45 °C, 2 min



Rate: 3 °C/min

Final Temp: 80 °C, 0 min

Detector: FID

Sensitivity: 128 x 10⁻¹² AFS

Injection Mode: Split



Components


1. Methanol
2. Propan-2-ol
3. Ethanol
4. sec-Butan-1-ol
5. n-Propanol + 2-Methyl-3-Buten-2-ol
6. d,l-3-Methyl-2-Butan-1-ol
7. Pentan-3-ol
8. d,l-2-Pentan-1-ol
9. n-Butanol
10. 2,4-Dimethyl, Pentan-3-ol
11. Hexan-3-ol
12. 2-Methyl Prop-2-en-1-ol
13. Crotyl Alcohol (2-Buten-1-ol)
14. Hexan-2-ol
15. Pentan-1-ol
16. 2-Methyl Pentan-1-ol
17. 2-Ethyl Butan-1-ol
18. Hexan-1-ol

ACI 03 | Analysis of 11 Organic Acids on BP20

| | | | |
|-------------------------|----------------------|-----------------|----------------------------|
| Column Part No.: | 054427 | | |
| Phase: | BP20, 0.25 µm film | Detector: | FID |
| Column: | 30 m x 0.25 mm ID | Sensitivity: | 32 x 10 ⁻¹² AFS |
| Initial Temp: | Isothermal at 155 °C | Injection Mode: | Split |



GC Columns and Applications



Components


1. 2,6-Xylenol
2. o-Cresol
3. Phenol
4. o-Ethylphenol
5. 2,5-Xylenol
6. p-Cresol
7. 2,4-Xylenol
8. m-Cresol
9. 2-iso Propylphenol
10. 2,3-Xylenol
11. 3,5-Xylenol + p-Ethylphenol

Notes: BP20 column completely resolves the three cresol isomers.

ACI 02 | Analysis of Organic Acids in Water on BP21

| | | | |
|-------------------------|-------------------|-----------------|----------------------------|
| Column Part No.: | 054477 | | |
| Phase: | BP21, 0.5 µm film | Final Temp: | 180 °C, 5 min |
| Column: | 30 m x 0.53 mm ID | Detector: | FID |
| Initial Temp: | 85 °C, 0 min | Sensitivity : | 64 x 10 ⁻¹² AFS |
| Rate: | 6 °C/min | Injection Mode: | On-Column |

Notes: On-column injection and the addition of a 0.03 M Oxalic acid (2%) to the injection solution increases the acidity of the column to allow lactic acid to be detected.




Components

1. Acetic Acid
2. Propanoic Acid
3. iso-Butyric Acid
4. n-Butyric Acid
5. iso-Valeric Acid
6. n-Valeric Acid
7. n-Caproic Acid
8. n-Heptanoic Acid
9. Lactic Acid






GC Columns and Applications



AMI 06 | Analysis of Aromatic Amines from Diazo Dyes on BPX35

| | |
|-------------------------|------------------------|
| Column Part No.: | 054701 |
| Phase: | BPX35 0.25 µm film |
| Azo Dyes standard: | 10 ppm solution in DCM |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp: | 50 °C, 2 min |
| Rate 1: | 15 °C to 240 °C |
| Rate 2: | 10 °C to 280 °C |
| Final Temp: | 280 °C, 25 min |
| Detector Type: | MS D |
| Carrier Gas: | He, 7.1 psi |
| Carrier Gas Flow: | 1.0 mL/min |

| | |
|--------------------------|----------------------------|
| Constant Flow: | On |
| Average Linear Velocity: | 36 cm/sec at 50 °C |
| Injection Mode: | Splitless |
| Purge on Time: | 1.0 min |
| Purge on (Split) | |
| Vent Flow: | 60 mL/min |
| Injection Volume: | 1 µL |
| Injection Temp: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part No.: | 092018 |




Components

- | | |
|---------------------------|---|
| 1. Indolin | 8. Unknown |
| 2. o-Toluidine | 9. 2-Naphthylamine |
| 3. 2,4-Diaminoanisole | 10. Unknown |
| 4. p-Chloroaniline | 11. 4-Aminodiphenyl |
| 5. p-Cresidine | 12. 4,4'-Oxydianiline |
| 6. 2,4,6-Trimethylaniline | 13. 4,4'-Diaminodiphenylmethane |
| 7. 4-Chlorotoluidine | 14. Benzidine |
| | 15. 3,3'-Dimethyl-4,4'-diaminodiphenylmethane |
| | 16. 3,3'-Dimethylbenzidine |
| | 17. 4,4'-Thiodianiline |
| | 18. 3,3'-Dichlorobenzidine |
| | 19. 4,4'-Methylene-bis(2-chloroaniline) |
| | 20. 3,3'-Dimethoxybenzidine |

AMI 03 | Analysis of Aromatic Amines on BP5

| | |
|-------------------------|-------------------|
| Column Part No.: | 054197 |
| Phase: | BP5, 1.0 µm film |
| Column: | 12 m x 0.53 mm ID |
| Initial Temp: | 60 °C, 0 min |
| Rate: | 10 °C/min |


| | |
|-----------------|-----------------|
| Final Temp: | 190 °C, 0 min |
| Detector: | FID |
| Sensitivity : | 128 x 10-12 AFS |
| Injection Mode: | Split |



Components

1. Pyridine
2. 2-Methyl Pyridine
3. Aniline
4. Benzylamine
5. o-Toluidine
6. m-Toluidine
7. 2,6-Dimethylaniline
8. 1,4-Phenyldiamine
9. Nicotine
10. Biphenylamine
11. Bibenzylamine

AMI 04 | Analysis of Amines on BP1



| | |
|-------------------------|-------------------|
| Column Part No.: | 054097 |
| Phase: | BP1, 3.0 µm film |
| Column: | 12 m x 0.53 mm ID |
| Initial Temp: | 70 °C |
| Rate: | 10 °C/min |
| Final Temp.: | 250 °C |
| Carrier Gas: | Nitrogen |
| Injection Volume: | 0.1 µL |

Components

1. Aniline
2. Decylamine
3. Dicyclohexylamine
4. Dodecylamine
5. Tetradecylamine




GC Columns and Applications

SOL 25 | Analysis of 22 Ketones on BPX35

| | | | |
|-------------------------|----------------------------|--------------------------|----------------------------|
| Column Part No.: | 054701 | Constant Flow: | On |
| Phase: | BPX35, 0.25 µm film | Average Linear Velocity: | 35 cm/sec at 40 °C |
| Sample: | 300 ppm in dichloromethane | Injection Mode: | Split |
| Column: | 30 m x 0.25 mm ID | Split Ratio: | 80:1 |
| Initial Temp: | 40 °C, 5 min. | Injection Volume: | 0.5 µL |
| Rate: | 10 °C/min to 170 °C | Injection Temp.: | 250 °C |
| Final Temp: | 170 °C, 5 min. | Liner Type: | 4 mm ID Single Taper Liner |
| Detector Type: | Mass Spectrometer | Liner Part Number: | 092017 |
| Carrier Gas: | He, 25.6 psi | Full Scan / SIM: | Full scan 45-450 |
| Carrier Gas Flow: | 1.6 mL/min. | | |


Components

- | | | |
|-------------------------|-------------------------|-------------------|
| 1. Acetone | 7. 3-Methyl-2-pentanone | 15. 3-Heptanone |
| 2. 2-Butanone | 8. 3-Hexanone | 16. 2-Heptanone |
| 3. 3-Methyl-2-butanone | 9. 2-Hexanone | 17. Cyclohexanone |
| 4. 2-Pentanone | 10. Mesityl oxide | 18. 2-Octanone |
| 5. 3-Pentanone | 11. 2-Methyl-3-hexanone | 19. 2-Nonanone |
| 6. 4-Methyl-2-pentanone | 12. Cyclopentanone | 20. 2-Decanone |
| | 13. 4-Methyl-2-hexanone | 21. 2-Undecanone |
| | 14. 5-Methyl-2-hexanone | 22. 2-Dodecanone |





GC Columns and Applications




SOL 33 | Analysis of Ketones on Thick Film BPX5

| | |
|-------------------------|----------------------------|
| Column Part No.: | 054123 |
| Phase: | BPX5, 1.0 µm film |
| Sample: | 300 ppm in dichloromethane |
| Column: | 60 m x 0.25 mm ID |
| Initial Temp.: | 40 °C, 5 min. |
| Rate 1: 1 | 0 °C/min to 80 °C |
| Rate 2: | 30 °C/min to 260 °C |
| Final Temp: | 260 °C, 4 min. |
| Detector Type: | FID |
| Detector Temp.: | 360 °C |

| | |
|--------------------------|----------------------------|
| Carrier Gas: | He, 27.6 psi |
| Carrier Gas Flow: | 1.9 mL/min. |
| Constant Flow: | On |
| Average Linear Velocity: | 35 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 100:1 |
| Injection Volume: | 0.4 µL |
| Injection Tem.: | 250 °C |
| Liner Type: | 4 mm ID Single Taper Liner |
| Liner Part Number: | 092017 |

Components

- | | | |
|-------------------------|-------------------------|-------------------|
| 1. Ethanol | 8. 3-Methyl-2-pentanone | 17. 2-Heptanone |
| 2. Acetone | 9. 3-Hexanone | 18. Cyclohexanone |
| 3. 2-Butanone | 10. 2-Hexanone | 19. 2-Octanone |
| 4. 3-Methyl-2-butanone | 11. Mesityl oxide | 20. 2-Nonanone |
| 5. 2-Pentanone | 12. Cyclopentanone | 21. 2-Decanone |
| 6. 3-pentanone | 13. 2-Methyl-3-hexanone | 22. 2-Undecanone |
| 7. 4-Methyl-2-pentanone | 14. 4-Methyl-2-hexanone | 23. 2-Dodecanone |
| | 15. 5-Methyl-2-hexanone | |
| | 16. 3-Heptanone | |





ALC 09 | Analysis of Phenols Mixture on BPX35




GC Columns and Applications

| | | | |
|------------------------|---------------------|---------------------------------|----------------------------|
| Column Part No: | 054701 | Constant Flow: | On |
| Phase: | BPX35, 0.25 µm film | Average Linear Velocity: | 35 cm/sec at 80 °C |
| Sample: | 200 ppm in methanol | Injection Mode: | Split |
| Column: | 30 m x 0.25 mm ID | Split Ratio: | 100:1 |
| Initial Temp: | 80 °C, 1 min | Injection Volume: | 1 µL |
| Rate 1: | 10 °C/min to 300 °C | Injection Temperature: | 250 °C |
| Final Temp: | 300 °C, 5 min | Liner Type: | 4 mm ID Single Taper Liner |
| Detector Type: | Mass Spectrometer | Liner Part No.: | 092017 |
| Carrier Gas: | He, 29.2 psi | Full Scan / SIM: | Full scan 45-450 |
| Carrier Gas Flow: | 1.7 mL/min. | | |





GC Columns and Applications




POL 06 | Analysis of Various Monomers on SolGel-WAX™

| | |
|-------------------------|--------------------------|
| Column Part No.: | 054796 |
| Phase: | SolGel-WAX, 0.25 µm film |
| Sample: | 250 ppm in Hexane |
| Column: | 30 m x 0.25 mm ID |
| Initial Temp: | 40 °C, 1 min. |
| Rate 1: | 10 °C/min to 250 °C |
| Final Temp: | 250 °C, |
| Detector Type: | FID |
| Detector Temp.: | 320 °C |
| Carrier Gas: | He, 16.6 psi |

| | |
|--------------------------|----------------------------|
| Carrier Gas Flow: | 1.6 mL/min. |
| Constant Flow: | On |
| Average Linear Velocity: | 35 cm/sec at 40 °C |
| Injection Mode: | Split |
| Split Ratio: | 80:1 |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Autosampler: | No |
| Liner Type: | 4 mm ID Single Taper Liner |
| Liner Part Number: | 092017 |

Components

- | | |
|--------------------------|---------------------------|
| 1. Ethyl acrylate | 6. Acrylic acid |
| 2. n-Butyl acrylate | 7. Hydroxypropyl acrylate |
| 3. n-Butyl methacrylate | 8. Methacrylic acid |
| 4. Styrene | 9. Unknown |
| 5. 2-Ethylhexyl acrylate | 10. Hydroxyethyl acrylate |
| | 11. Lauryl Methacrylate |



POL 01 | Analysis of Unreacted Monomers in Latex on BP20


| | |
|-------------------------|-------------------|
| Column Part No.: | 054488 |
| Phase: | BP20, 1.0 µm |
| Column: | 25 m x 0.53 mm ID |
| Initial Temp.: | 40 °C, 2 min |
| Rate: | 10 °C/min |

| | |
|-----------------|-----------------|
| Final Temp.: | 230 °C, 5 min |
| Injector Cond.: | Split, 280 °C |
| Detector: | FID, 280 °C |
| Carrier Gas: | Hydrogen, 4 psi |

Note: This was performed by heated headspace analysis.

Components

- Vinyl Acetate
- Ethyl Acrylate
- Monomethyl Methacrylate
- Butyl Acrylate
- Butyl Methacrylate
- Styrene
- Di-methylamino Ethyl-methacrylate
- 2-Ethyl Hexylacrylate
- Octanol
- Unknown
- 2-(acetoacetoxy) Ethyl Methacrylate
- Dibutyl Maleate
- Dicyclopentenloxyethyl Methacrylate





POL 05 | Pyrolysis of Styrene-isoprene Copolymer Pyrolysis of Polystyrene on BP1

| | |
|-------------------------|-----------------------|
| Column Part No.: | 054053 |
| Phase: | BP1, 1.0 µm |
| Column: | 25 m x 0.22 mm ID |
| Initial Temp.: | 40 °C, 1 min |
| Rate: | 10 °C/min |
| Final Temp.: | 140 °C |
| Detector: | FID |
| Pyrolysis Temp.: | 550 °C |
| Carrier Gas: | H ₂ 10 psi |

| | |
|-------------------------|------------------------|
| Column Part No.: | 054065 |
| Phase: | BP1, 0.5 µm |
| Column: | 25 m x 0.32 mm ID |
| Initial Temp.: | 40 °C, 1 min |
| Rate: | 10 °C/min |
| Final Temp.: | 130 °C |
| Detector: | FID |
| Pyrolysis Temp.: | 800 °C |
| Carrier Gas: | H ₂ , 5 psi |




GC Columns and Applications



SOL 21 | Analysis of a Common Solvent Mixture on BP624


| | |
|-------------------------|--------------------------------|
| Column Part No.: | 054832 |
| Phase: | BP624, 1.8 µm film |
| Alcohol mix: | 1000 ppm in Dimethyl Sulfoxide |
| Column: | 30 m x 0.32 mm ID |
| Initial Temp: | 32 °C, 9 min. |
| Rate: | 30 °C/min to 190 °C |
| Final Temp: | 190 °C, 0 min. |
| Detector Type: | FID |
| Carrier Gas: | He, 9.6 psi |
| Carrier Gas Flow: | 2.2 mL/min. |

| | |
|--------------------------|----------------------------|
| Constant Flow: | On |
| Average Linear Velocity: | 34 cm/sec at 32 °C |
| Injection Mode: | Split |
| Split Ratio: | 100:1 |
| Injection Volume: | 0.2 µL |
| Injection Temperature: | 250 °C |
| Autosampler: | No |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |






GC Columns and Applications



ENV 52 | Industrial Solvents on SolGel-WAX™

| | | | |
|--------------------------|--|----------------------------|--------------------|
| Column Part No. | 054797 | Constant Flow: | On |
| Phase: | SolGel-WAX™, 0.5 µm film 30 m x 0.32 mm ID | Pressure: | 8.4 psi |
| Split / Splitless | | Column Flow: | 1.84 mL/min |
| Injector Temp: | 240 °C | Linear Velocity: | 30 cm/sec at 35 °C |
| Injection Volume: | 0.1 µL | Initial Temp.: | 35 °C |
| Autosampler Syringe: | 0.5 µL Removable Needle Part No. 000410 | Initial Time: | 3 min |
| Septa: | Auto-Sep T™ Part No. 041882 | Rate 1: | 15 °C/min |
| Injection Type: | Split | Final Temp. 1: | 230 °C |
| Purge On Time: | NA | Hold Time: | 4 min |
| Purge On (Spilt) Vent: | 150 mL/min | Run Time: | 20.00 min |
| Split Ratio: | 83 to 1 | Detector Parameters | |
| Liner Type: | Single taper Part No. 092017 | Detector Type: | FID at 270 °C |
| Carrier Gas: | He | | |

Sample Description: Industrial solvents mix, 25 to 50 ng per component on column




Components

1. Acetone
2. Ethyl acetate
3. Methyl ethyl ketone
4. Contaminant
5. iso-Propanol
6. Ethanol
7. Methyl isobutyl ketone
8. Toluene
9. Butyl acetate
10. iso-butanol
11. Propylene glycol monomethyl ether
12. n-Butanol
13. Ethyl benzene
14. p-Xylene
15. m-Xylene
16. o-Xylene
17. Butyl Cellosolve acetate
18. Cyclohexanone
19. Butyl Cellosolve
20. Butyl glycol acetate
21. Hexyl Cellosolve
22. Isophorone
23. Butyl Carbitol
24. Benzyl alcohol

SOL 04 | Acrylic Acid/Acrylamide Analysis on BP21

| | |
|-------------------------|-------------------|
| Column Part No.: | 054473 |
| Phase: | BP21, 0.5 µm film |
| Column: | 12 m x 0.53 mm ID |
| Initial Temp: | 75 °C, 0.5 min |
| Rate: | 10 °C/min |
| Final Temp: | 150 °C |
| Detector: | FID, 280 °C |
| Injection Mode: | On-Column |
| Carrier Gas: | He, 6 psi |


Notes: When response of acrylic acid is low, removal of 30 cm from the front of the column will correct this loss. On-column injection is recommended or polymerization of acrylic acid may occur.



Components


1. Acrylic acid (25 ng)
2. Acrylamide (10 ng)

TP-0138-C | Analysis Of Polybrominated Diphenyl Ethers on BP1



SGE would like to thank the Japan Food Research Centre for evaluating the BP1 column, SGE Japan and Chemicals Evaluation and Research Institute, Japan Toshiyuki KATAOKA, Masahiro AKIBA and Shinnichi KUDO.

TP-0138-C | Analysis Of Polybrominated Diphenyl Ethers on BPX5



SGE would like to thank SGE Japan and Chemicals Evaluation and Research Institute, Japan Toshiyuki KATAOKA, Masahiro AKIBA and Shinnichi KUDO.




GC Columns and Applications







GC Columns and Applications



TP-0138-C | Analysis Of A Mixture Of Pentachlorobiphenyls on HT8-PCB



The separation of a mixture of pentachlorobiphenyls using an HT8-PCB column. Elution order calculated for the 5CBs from structure activity relationships based on coplanarity and confirmation, steric factors and electron density show a high correlation with experimental results.

SGE would like to thank T. Nakano, C. Matsumura and M. Tsurukawa at Hyogo Prefectural Institute of Public Health and Environmental Sciences, for providing the PCBs on HT8-PCB data.

TP-0138-C | Analysis Of A Mixture Of PBDD, PCDD And PBDF on BPX70





Figure 4. The separation of a mixture of PBDD and PBDF on a BPX70 column. The mixture was separated using the π - π interaction between the compounds and the cyano phase of the BPX70 column.


SGE would like to thank Toshiyuki Kataoka, Masahiro Akiba and Shinnichi Kudo of the Chemicals Evaluation and Research Institute, Japan, and SGE Japan, for providing the chromatograms of PBDEs on the ENV-5 and BPX70 columns.

PHA 14 | Analysis of Drugs of Abuse on BPX35

| | | | |
|-------------------------|---------------------|--------------|---------------|
| Column Part No.: | 054711 | Temp 2: | 200 °C |
| Phase: | BPX35, 0.25 µm film | Rate 2: | 7 °C/min |
| Column: | 25 m x 0.22 mm ID | Temp 3: | 295 °C |
| Initial Temp.: | 80 °C | Rate 3: | 20 °C/min |
| Rate 1: | 15 °C/min | Final Temp.: | 340 °C, 6 min |



GC Columns and Applications




PHA 09 | Analysis of Tricyclic Antidepressants on BPX35

| | |
|-------------------------|-------------------|
| Column Part No.: | 054711 |
| Phase: | BPX35, 0.25 µm |
| Column: | 25 m x 0.22 mm ID |
| Initial Temp.: | 210 °C, 1 min |
| Rate: | 5 °C/min |
| Final Temp.: | 280 °C |
| Carrier Gas: | Helium, 150 kpa |
| Injection Mode: | Split (20:1) |
| Detector: | FID, 380 °C |

Components

1. Amitriptyline
2. Trimipramine
3. Nortriptyline
4. Doxepin
5. Desipramine




Note: BPX35 is a low bleed, chemically inert phase which allows trace analysis to occur.






GC Columns and Applications



PHA 19 | Analysis of a Variety of Antidepressant and Anticonvulsant Drugs on BPX50

| | |
|-------------------------|----------------------|
| Column Part No.: | 054751 |
| Phase: | BPX50, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Sample: | 5-10 ppm in methanol |
| Initial Temp: | 150 °C, 0.5 min |
| Rate 1: | 10 °C/min to 180 °C |
| Rate 2: | 1.5 °C/min to 220 °C |
| Rate 2: | 30 °C/min to 260 °C |
| Final Temp: | 260 °C, 5 min |
| Detector Type: | FID |
| Detector Temp.: | 320 °C |
| Carrier Gas: | He, 25.7 psi |

| | |
|-----------------------------|----------------------------|
| Carrier Gas Flow: | 1.8 mL/min. |
| Constant Flow: | On |
| Average Linear Velocity: | 35 cm/sec at 40 °C |
| Injection Mode: | Splitless |
| Purge on Time: | 0.5 min |
| Purge on (Split) Vent Flow: | 60 mL/min |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Liner Type: | 4 mm ID Single Taper Liner |
| Liner Part Number: | 092017 |
| Full Scan / SIM: | Full scan 45-450 |



Components

1. Benzothiazole
2. Selegeline
3. Methsuximide
4. Ketamine
5. Unknown
6. Brompheniramine
7. Imipramine
8. Mianserin
9. Dothiepin
10. Carbamazepine
11. Prazepam
12. Diltiazem
13. Pholcodine


PHA 13 | Analysis of Blood Alcohol on BP20

| | |
|-------------------------|----------------------------|
| Column Part No.: | 054442 |
| Phase: | BP20, 1.0 µm film |
| Column: | 25 m x 0.32 mm ID |
| Initial Temp: | Isothermal, 60 °C |
| Detector: | FID |
| Sensitivity: | 64 x 10 ⁻¹² AFS |
| Injection Mode: | Split |

Note: The BP20 column allows the use of aqueous solutions.

Components


1. Acetone
2. Ethyl Acetate
3. Methanol
4. iso-Propanol
5. Ethanol
6. n-Propanol



PHA 06 | Analysis of Basic Drug Screen on BPX5 (10-20 ng/component)

| | | | |
|-------------------------|---------------------|--------------|------------------------|
| Column Part No.: | 054131 | | |
| Phase: | BPX5, 1.0 µm | Final Temp.: | 310 °C |
| Column: | 25 m x 0.53 mm I.D. | Detector: | FID |
| Initial Temp.: | 120 °C | Injector: | Split, 240 °C |
| Rate: | 10 °C/min | Carrier Gas: | H ₂ , 2 psi |

Note: The low bleed nature of the BPX5 allows trace analysis to be performed.




Components

1. Methamphetamine
2. Phendimetrazine
3. Phencyclidine
4. Mepivocaine
5. Methaqualone
6. Amitriptyline
7. Codeine
8. Diazepam
9. Fentanyl



GC Columns and Applications

PHA 08 | Underivatized Steroid Analysis on BPX5



| | |
|-------------------------|----------------------------|
| Column Part No.: | 054113 |
| Phase: | BPX5, 0.25 µm |
| Column: | 25 m x 0.22 mm ID |
| Initial Temp.: | 180 °C |
| Rate: | 8 °C/min |
| Final Temp.: | 350 °C, 10 min |
| Detector: | FID |
| Sensitivity: | 32 x 10 ⁻¹² AFS |
| Injection Mode: | Split |
| Carrier Gas: | H ₂ , 10 psi |

Note: The low bleed nature of the BPX5 allows trace analysis to be performed.


Components

1. Testosterone
2. Pregnenolone
3. Progesterone
4. Cholesterol






GC Columns and Applications



PHA 15 | Analysis of Horse Racing Test Mix on BPX5

| | |
|-------------------------|---------------------|
| Column Part No.: | 054101 |
| Phase: | BPX5, 0.25 µm film |
| Column: | 30 m x 0.25 mm ID |
| Horse Racing standard*: | 10 ppm in methanol |
| Initial Temp: | 75 °C, 2 min |
| Rate 1: | 15 °C/min to 300 °C |
| Rate 2: | 20 °C/min to 320 °C |
| Final Temp: | 320 °C, 8 min. |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 14.5 psi |
| Carrier Gas Flow: | 1.5 mL/min |


| | |
|--------------------------|----------------------------|
| Constant Flow: | On |
| Average Linear Velocity: | 45 cm/sec at 75 °C |
| Injection Mode: | Splitless |
| Purge on Time: | 0.5 min |
| Purge on (Split) | |
| Vent Flow: | 60 mL/min |
| Injection Volume: | 1 µL |
| Injection Temperature: | 250 °C |
| Liner Type: | 4 mm ID Double Taper Liner |
| Liner Part Number: | 092018 |



Components

- | | |
|--------------------|-------------------|
| 1. Metronidazole | 6. Trimipramine |
| 2. Amylobarbitone | 7. Phenytoin |
| 3. Pentobarbitone | 8. Diazepam |
| 4. Caffeine | 9. Nordiazepam |
| 5. Diphenhydramine | 10. Diphenoxylate |

PHA 03 | Analysis of Alkaloids on BP5



| | |
|-------------------------|-----------------------------|
| Column Part No.: | 054198 |
| Phase: | BP5, 1.0 µm film |
| Column: | 25 m x 0.53 mm ID |
| Initial Temp.: | 200 °C, 0 min |
| Rate: | 25 °C/min |
| Final Temp: | 300 °C, 0 min |
| Detector: | FID |
| Sensitivity: | 128 x 10 ⁻¹² AFS |
| Injection Mode: | Split |

Note: A 0.53 mm ID column can be used to screen samples rapidly.


Components

1. Cocaine
2. Codeine
3. Morphine
4. Quinine

PHA 10 | Underivatized Barbiturates on BP5



GC Columns and Applications




| | |
|-------------------------|---------------------|
| Column Part No.: | 054197 |
| Phase: | BP5, 1.0 µm |
| Column: | 12 m x 0.53 mm I.D. |
| Temp: | 195 °C |
| Carrier Gas: | Hydrogen |
| Carrier Flow: | 10 mL/min |
| Injection Volume: | 0.1 µL |

Note: A 0.53 mm ID column can be used to screen samples rapidly.

Components

1. Barbital
2. Butabarbital
3. Amobarbital
4. Pentabarbital
5. Secobarbital
6. Hexabarbital

PHA 04 | Analysis of Sedatives/Hypnotics on BP1



| | |
|-------------------------|------------------------------|
| Column Part No.: | 054087 |
| Phase: | BP1, 1.0 µm film |
| Column: | 25 m x 0.53 mm ID |
| Initial Temp.: | 180 °C, 0 min |
| Rate: | 10 °C/min |
| Final Temp.: | 250 °C, 3 min |
| Detector: | FID |
| Sensitivity: | 1024 x 10 ⁻¹² AFS |
| Injection Mode: | Split |

Components

1. Allobarbital
2. Aprobarbital
3. Diphenhydramine
4. Mephobarbital
5. Methapyrilene
6. Chlorpheniramine







GC Columns and Applications

PHA 01 | Analysis of Acid/Neutral Drugs on BPX35

| | | | |
|-------------------------|-------------------|-----------------|---------------|
| Column Part No.: | 054711 | | |
| Phase: | BPX35, 0.25 µm | Final Temp.: | 300 °C, 5 min |
| Column: | 25 m x 0.22 mm ID | Carrier Gas: | He, 150 kpa |
| Initial Temp.: | 100 °C, 1 min | Injection Mode: | Split, (20:1) |
| Rate: | 10 °C/min | Detector: | FID, 380 °C |

Note: BPX35 is a low bleed column with a maximum temperature of 360 °C. Very compatible with GC/MS systems.



Components

- (Concentration: 200-400 ug/ml)
1. Ethosuximide
 2. Barbitol
 3. Aprobarbital
 4. Butabarbital
 5. Amobarbital
 6. Pentobarbital
 7. Secobarbital
 8. Meprobamate
 9. Carisoprodal
 10. Glutethimide
 11. Phenobarbital
 12. Methaqualone
 13. Primidone

