

Manual Crimpers and Decappers

Hand Operated Vial Crimper

Use to crimp vials and bottles in one quick, easy motion. Built-in adjustable stop allows the user to repeatedly achieve an optimum seal. Solid construction, combined with quality components manufactured to precise tolerances, provide consistent, dependable performance. Available for use with a wide variety of vial closure styles and sizes ranging from 8mm to 20mm.



Description	Part Number
8mm Hand Operated Crimper	9300-08
11mm Hand Operated Crimper	9300-11
13mm Hand Operated Crimper	9300-13
20mm Hand Operated Crimper	9300-20

Hand Operated Plier Type Decapper

Plier type decappers are a convenient, low cost hand tool which allows crimped aluminium seals to be removed quickly and easily from vials and bottles. Simple operation- Just squeeze the cap and twist off.

Four standard sizes to choose from

Sizes: 8mm, 11mm, 13mm, 20mm



Description	Part Number
8mm Plier Type Decapper	9310-08
11mm Plier Type Decapper	9310-11
13mm Plier Type Decapper	9310-13
20mm Plier Type Decapper	9310-20

Hand Operated Plier Type Decapper

Hand operated crimp type decappers are used to remove the crimp seal from vials and bottles in one smooth motion. Sizes range from 8mm to 20mm. Solid construction combined with quality materials and components provide a sturdy and reliable hand tool.



Description	Part Number
8mm Hand Operated Crimp Type Decapper	9320-08
11mm Hand Operated Crimp Type Decapper	9320-11
13mm Hand Operated Crimp Type Decapper	9320-13
20mm Hand Operated Crimp Type Decapper	9320-20